
Combur¹⁰ Test

cobas
V 100

REF 04510062 171

REF 04510089 056

English

Intended use: Ten-patch test strip for the semi-quantitative determination of specific gravity, pH, leukocytes, nitrite, protein, glucose, ketone bodies, urobilinogen, bilirubin, and blood in urine by visual reading.

For professional use only.

Summary:

Urine test strips are used to measure certain constituents in urine which are significant of renal, urinary, hepatic and metabolic disorders.

Test principle

Specific gravity (SG): The test detects the ion concentration of the urine. In the presence of cations, protons are released by a complexing agent and produce a color change in the indicator bromothymol blue from blue via blue-green to yellow.

pH: The test paper contains the indicators methyl red, phenolphthalein and bromothymol blue and reacts specifically with H⁺-ions. The most frequent pH values of fresh urine from healthy subjects lie between 5 and 6.

Leukocytes (LEU): The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye. Bacteria, trichomonads or erythrocytes present in the urine do not affect the reaction.

Nitrite (NIT): The test is based on the principle of the Griess test and is specific for nitrite. The reaction reveals the presence of nitrite and hence indirectly nitrite-forming bacteria in the urine by a pink-to-red-coloration of the test patch. Even a slight pink coloration is indicative of significant bacteruria.

Protein (PRO): The test is based on the principle of the protein error of a pH indicator. It is particularly sensitive to albumin. An elevated pH (up to 9) does not affect the test.

Glucose (GLU): The glucose determination is based on the specific glucose-oxidase/peroxidase reaction (GOD/POD method). The test is independent of the pH and specific gravity of the urine and is not affected by the presence of ketone bodies.

Ketones (KET): This test is based on the principle of Legal's test and is more sensitive to acetoacetic acid than to acetone.

Urobilinogen (UBG): A stable diazonium salt reacts almost immediately with urobilinogen to give a red azo dye. The test is specific for urobilinogen and is not susceptible to the interfering factors known to affect the Ehrlich's test.

Bilirubin (BIL): The test is based on the coupling of bilirubin with a diazonium salt. Even the slightest pink coloration constitutes a positive, i.e. pathologic, result. Other urinary constituents produce a more or less intense yellow coloration.

Blood (ERY/Hb): The peroxidase-like action of hemoglobin and myoglobin specifically catalyzes the oxidation of the indicator by means of the organic hydroperoxide contained in the test paper to give a blue-green coloration.

Reagents: Each test contains per 1 cm² test patch area the following:

Specific gravity: Ethyleneglycol-bis(diaminoether)tetraacetic acid 182.8 µg; bromothymol blue 36 µg

pH: Bromothymol blue 13.9 µg; methyl red 1.2 µg; phenolphthalein 8.6 µg

Leukocytes: Indoxylcarboxylic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg

Nitrite: 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg

Protein: 3',3'',5'-tetrachlorophenol-3,4,5,6-tetrabromosulfophthalein 13.9 µg

Glucose: 3',3'',5,5'-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U

Ketones: Sodium nitroprusside 157.2 µg

Urobilinogen: 4-methoxybenzene-diazonium-tetrafluoroborate 67.7 µg

Bilirubin: 2,6-dichlorobenzene-diazonium-tetrafluoroborate 16.7 µg

Blood: 3',3'',5,5'-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

Precautions and warnings:

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

The stopper of the test strip vial contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

Reagent handling:

Test strips are ready for use.

Storage and stability:

Store the package at 2-30 °C. The test strips are stable up to the expiration date specified on the box, when stored in the original container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation: Use only a fresh midstream urine sample¹ that has not been centrifuged.

For specimen collection and preparation only use suitable tubes or collection containers.

Use only clean, well-rinsed vessels to collect urine.

False-positive readings for blood and glucose can result from residues of strongly oxidizing disinfectants in the specimen collection vessel. The urine specimen should not stand for more than 2 hours before testing. In case of longer standing, store the specimen refrigerated at 4°C.

At the time of testing, ensure the sample is at room temperature and mix the sample thoroughly before use. Longer standing times can lead to false results owing to the following influences: Proliferation of bacteria; A rise in pH due to ammonia formed as a result of bacterial degradation of urea.

Do not add preservatives to the urine.

Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these two parameters². Drugs that turn red in an acid environment (e.g. phenazopyridine) may produce false-positive readings or reddish colorations on the test patches for nitrite, protein, urobilinogen and bilirubin.

Materials provided

- REF 04510062171, package with 100 test strips

Materials required (but not provided)

- General laboratory equipment

Assay:

For optimum performance of the assay follow the directions given in this document.

- Use fresh urine that has not been centrifuged. Thoroughly mix the urine sample. The sample should be at room temperature when the test is performed and should not have been standing for more than 2 hours.
- Take a test strip out of the container. Close the container again with the original desiccator stopper immediately after removal of the strip. This is important as otherwise the test areas may become discolored due to environmental influences such as moisture or nitrite gases in the air. Incorrect results may be obtained.
- Briefly (about 1 second) dip the test strip into the urine making sure that all test areas are moistened.
- When withdrawing the test strip, wipe the edge against the rim of the vessel to remove excess urine.
- After 60 seconds (60-120 seconds for the leukocyte test area), compare the reaction colors of the test areas with the colors on the label and assign always the value of the nearest color block.
- Compare the 10th (blood) test area with both color scales as separate color scales are given for erythrocytes and hemoglobin.

Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

Quality control:

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

For quality control, use commercially available urine controls, or other suitable control material. Important note for reporting results (for professional users).

According to the regulations from the German Medical Association for quality assurance of medical laboratory analyses dated 11/23/2007, the decision to classify a laboratory test result to either B1 or B2 depends on the way the test results are expressed in the report (scale level).

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