

ALPHA AMYLASE SL

Kinetic Method (CNP3)
Liquid Reagent ready to use

REF. 5501 6x 10 ml
REF. 5502 3x 50 ml
REF. 5503 3x 10 ml



Azienda certificata DNV



INTENDED USE

Quantitative determination of α -Amylase in serum and plasma.

PRINCIPLE

The α -Amylase hydrolyzes 2-chloro-4-nitrophenyl- α -D-maltotriose (CNP3) into 2-chloro-4-nitrophenyl- α -D-maltoside (CNP2), maltotriose (G3), glucose and 2-chloronitrophenol. The absorbance change in unit time measured at 405 nm is proportional to the enzyme activity in the sample.

SAMPLE

Serum, heparinized plasma, urine.

Do not use other anticoagulants like EDTA, citrate and oxalate as they inhibit the enzyme. Avoid hemolyzed samples.

Dilute urine 1:3 with saline solution.

Amylase activity is stable one week at 15-25°C and some months at 2-8°C.

KIT COMPONENTS

Reagent (A) AMY Volume = 10/50 ml	Good buffer pH 6	100 mmol/l
	CNP3	3.1 mmol/l
	Sodium chloride	10 mmol/l
	Calcium acetate	1 mmol/l
	Sodium azide	15 mmol/l

The reagent is stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not freeze. Once opened reagent is stable for 2 months at 2-8°C if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent, bring to room temperature (15-25°C) before use.

PRECAUTIONS AND WARNINGS

In accordance with the regulations, the reagent is classified: **Xn – Harmful**.

R22 – Harmful if swallowed **R52/53** – Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

S60 – This material and its container must be disposed of as hazardous waste.

S61 – Avoid release to the environment.

Handle with care, avoiding skin contact and swallow.

Use the normal precautions required in the laboratory.

PROCEDURE

Wavelength:	405 nm
Lightpath :	1 cm
Temperature:	37°C
Reading:	against distilled water
Method:	Increasing Kinetic
Sample/Reagent:	1/50

pipette:

Reagent (A)	1000 μ l
sample	20 μ l

Mix, incubate at 37°C for 1 minute, read the initial absorbance against water. Perform 3 readings at 60 seconds intervals. Calculate the average value of the absorbance variations per minute (ΔA /min).

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

Perform calculation in units per litre, multiplying the ΔA /min by the factor:
activity in U/L: ΔA /min \times 3517

EXPECTED VALUES

Serum/plasma: 25 – 110 U/L

Urine: < 480 U/L

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 1 U/L

Linearity: the method is linear up to 2350 U/L. For higher values, dilute the sample 1:2 and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (U/l)	42.90	184.2	497.5
DS	0.64	2.25	5.64
CV %	1.49	1.22	1.13

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (U/l)	44.15	181.2	524.5
DS	0.48	2.82	2.22
CV %	1.09	1.56	0.42

Interferences: bilirubin does not interfere up to 30 mg/dl. Triglycerides up to 1500 mg/dl do not interfere. Hemoglobin up to 500 mg/dl does not interfere. Glucose up to 500 mg/dl does not interfere.

Correlation against a reference method: $Y = 0.4582x + 1.089$ $r = 0.999$

REFERENCES

1. Lorentz, K., Clin. Chem. Clin. Biochem. 17, 499 (1979).
2. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
3. Young, D.S., et al., Clin. Chem. 21:1D (1975).