

# CK-MB

Kinetic UV Method.

REF. 0083 40x 10 ml



DNV CERTIFIED COMPANY  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012



## INTENDED USE

Quantitative determination of CK-MB in serum.

## PRINCIPLE

CK-MB consists of the subunits CK-M and CK-B. Specific antibodies against CK-M inhibits completely CK-MM activity (main part of the total CK activity) and CK-M subunit of CK-MB. Therefore only CK-B activity is measured, which is half of the CK-MB activity.

The activity of CK-MB is obtained multiplying the activity of CK-B by 2

## SAMPLE

Serum, plasma with heparin. Do not use hemolyzed samples.

CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes. The CK is subject to thermal denaturation, rapidly cool the sample at 4°C after collection.

## KIT COMPONENTS

Reagent (A) CK-MB Volume = 40 ml	Imidazole buffer, pH 6.7	150 mmol/l
	D-Glucose	25 mmol/l
	N-Acetyl-L-cysteine	25 mmol/l
	Magnesium acetate	12.5 mmol/l
	NADP	2.52 mmol/l
	EDTA	2.02 mmol/l
Reagent (B) CK-MB Volume = 10 ml	Hexokinase	≥6800 U/l
	Anti-CK-M antibody with inhibitory capacity	2000 U/l
	Creatine phosphate	250 mmol/l
	ADP	15.2 mmol/l
	AMP	25 mmol/l
	Diadenosine pentaphosphate	10 mmol/l
	G-6-PDH	≥8800 U/l

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not use over expiry data.

Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided. Keep bottles closed when not in use.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

## REAGENT PREPARATION

Add 1 Volume of Reagent (B) to 4 Volumes of Reagent (A). Ratio 4+1.

The working solution (A+B) is stable 1 week at 2-8°C.

## PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## PROCEDURE

Wavelength:	340 nm
Lightpath :	1 cm
Temperature:	37°C
Reading:	against distilled water
Method:	increasing kinetic

pipette:

Working solution (A+B)	1000 µl
sample	50 µl

Mix, incubate at 37°C for 10 minutes, read the initial absorbance(A) of the sample , start the stopwatch and read again after 5 minutes (A<sub>2</sub>). Calculate the difference between absorbances ΔA=A<sub>2</sub> - A<sub>1</sub> .

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

## RESULTS CALCULATION

Perform calculation in Units per litre, multiplying the ΔA by the factor as it is indicated:

$$\text{CK-MB Activity (U/L):} \quad \Delta A \times 4500 (*)$$

Percentage of CK-MB in the sample:

$$\% \text{ CK-MB} = \text{CK-MB} / \text{Total CK} \times 100$$

to calculate **Total CK activity** use the kits CK-NAC SL (Ref. 0018-4054-4052)

(\* Factor calculated in our laboratories. We recommend the use of CK-MB Calibrator (Ref. 6021 - 1x2 ml) to verify that this factor is correct for your test system.

## EXPECTED VALUES

CK-MB	< 25 U/l
% CK-MB	< 6 %

The probability of myocardial infarction is high under the following conditions:

### Total CK

Men	> 171 U/l
Women	> 145 U/l

CK-MB	> 25 U/l
% CK-MB	6 – 25 %

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: CK-MB Control (Ref. 6020).

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is: 2 U/L.

**Linearity:** the method is linear up to 600 U/L. For higher values, dilute the sample 1:10 and multiply the result by 10.

### Precision intra-assay:

	Level 1	Level 2
Mean (U/L)	25	66
CV %	10.3	4.5

### Precision inter-assay:

	Level 1	Level 2
Mean (U/L)	25	74
CV %	9.8	2.6

**Interferences:** No significant interference up to 8 mmol/l of triglycerides. Glucose does not interfere up to 7g/l. Hemoglobin does not interfere up to 6g/l .

**Correlation against a reference method:** Y = 1.018x+0.308 r = 0.99

## REFERENCES

1. Abbot B et al. Creatinine Kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116
2. Gerhardt W et al. Creatine Kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young, D.S., Effects of drugs on Clinical Lab. Tests, 4<sup>th</sup> Ed., AACC Press, 1995.
4. Young, D.S., Effects of disease on Clinical Lab. Tests, 4<sup>th</sup> Ed., AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd Ed. AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd Ed. AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la creatinine kinase dans la serum humain. Ann. Biol. Clin., 40, (1482),87.

## Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 05111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) – web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

008207

Ed. 2015/12 rev. 00