

CK-NAC SL

Kinetic Method UV

Liquid reagents ready to use

REF. 0018 5x 10 ml
REF. 4054 2x 50 ml
REF. 4052 2x100 ml



Azienda certificata DINV



INTENDED USE

Quantitative determination of Creatine kinase, CK, in serum and plasma.

PRINCIPLE

Creatin kinase catalyzes the hydrolysis reversible reaction of creatin phosphate in creatinine and ATP. ATP, in presence of exokinase and D-glucose, is transformed into ADP and D-glucose-6-phosphate, which reduces NADP⁺ in NADPH, H⁺ and D-gluconate-6-phosphate in presence of glucose-6P dehydrogenase. The increase absorbance at 340 nm in time unit, due to NADPH formation, defines the CK activity in the sample.

SAMPLE

Serum, plasma with heparin. Do not use hemolyzed samples.

A slight degree of hemolysis, up to 200 mg/dl of Hb, can be tolerated.

CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes.

KIT COMPONENTS

Component	Quantity
Reagent (A) CK Volume = 10/40/80 ml	Good buffer 125 mmol/l Magnesium acetate 12 mmol/l EDTA 2 mmol/l D-glucose 25 mmol/l N-acetyl-L-cysteine 25 mmol/l NADP 2.5 mmol/l HK-hexokinase 6500 U/l
Reagent (B) CK Volume = 10/20 ml	ADP 15 mmol/l AMP 25 mmol/l Diadenosine-5-phosphate 103 mmol/l G-6-PDH 8800 U/l Creatine phosphate 250 mmol/l

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not freeze. Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent, bring to room temperature (15-25°C) before use.

For use as monoreagent: add a part of Reagent (B) to 4 parts of Reagent (A).

The working solution (A+B) is stable 2 days at 15-25°C and 2 weeks at 2-8°C.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength:	340 nm (334 – 365)
Lightpath :	1 cm
Temperature:	37°C
Reading:	against distilled water
Method:	Increasing Kinetic
Sample/Reagent:	1/25

Use as monoreagent:

pipette:

Working solution (A+B) 1000 µl

sample 40 µl

Mix, incubate at 37°C for 1 minute, read the initial absorbance against water. Perform 3 readings at 60 seconds intervals. Calculate the average value of the absorbance variations per minute. (ΔA/min).

Use as bireagent:

pipette:

Reagent (A) 1000 µl

sample 50 µl

mix and after 1 minute add:

Reagent (B) 250 µl

Mix, incubate at 37°C for 1 minute, read the initial absorbance against water. Perform 3 readings at 60 seconds intervals.

Calculate the average value of the absorbance variations per minute. (ΔA/min).

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

Perform calculation in Units per litre, multiplying the ΔA/min by the factor as it is indicated:

Activity in U/L: ΔA/min x 4127

EXPECTED VALUES

Men: 24 - 204 U/L

Women: 24 - 173 U/L

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 1 U/L

Linearity: the method is linear up to: 1200 U/L. For higher values, dilute the sample 1:10 and multiply the result by 10.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (U/l)	73.41	156.9	499.1
DS	0.590	0.738	0.876
CV %	0.80	0.47	0.18

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (U/l)	78.93	157.3	502.4
DS	0.523	1.059	1.430
CV %	0.66	0.67	0.29

Interferences: do not use anticoagulants different from heparin because they inhibit CK activity. Bilirubin does not interfere up to 5 mg/dl. Hemoglobin does not interfere up to 0.2 g/dl.

Correlation against a reference method: Y = 0.9712x - 2.5338 r = 0.9996

REFERENCES

1. Mathieu M. et al., Ann. Biol.Clin., 40, 99 (1982).
2. Young, D.S., et al., Clin. Chem. 21:1D (1975).
3. Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby Ed. 1989.
4. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).