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Kit Controllo HbA_{1c} – REF. 6744.

INTENDED USE

This product is for the calibration in the quantitative determination of Glycated Hemoglobin (HbA_{1c}) in blood.

PRODUCTS CHARACTERISTICS

The lyophilized HbA_{1c} calibrators are a hemolysate prepared from human erythrocytes.

KIT COMPONENTS

HbA _{1c} Calibrator 1 Volume = 0.5 ml	Human blood. Contains preservatives. The HbA_{1c} concentration is indicated on the label
HbA _{1c} Calibrator 2 Volume = 0.5 ml	
HbA _{1c} Calibrator 3 Volume = 0.5 ml	
HbA _{1c} Calibrator 4 Volume = 0.5 ml	

PREPARATION

Reconstitute each calibrator with 0.5 ml of deionized water.
Gently mix for 10 minutes, or until all material has dissolved.

PROCEDURE

The lyophilized HbA_{1c} calibrators should be treated in the same manner as specimens, hemolysate procedure included.

The calibrators produce a calibration curve that will be stable for at least 7 days on most analyzers.

STORAGE AND STABILITY

Store at 2-8°C. The calibrators are stable until expiration date if sealed tightly, protected from light and heat.

The reconstituted calibrators are stable 30 days at 2-8°C and 3 months at -20°C. Do not freeze and thaw more than once.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HbsAg, HCV and antibody to HIV 1/2. However handle cautiously as potentially infectious.

MATERIAL REQUIRED AND NOT PROVIDED

Kit Direc Glycate Hemoglobin – REF. 6738.

Kit HbA_{1c} Control – REF. 6744.