# ISETROL

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**Intended use**

ISETROL, an electrolyte control is intended for use as control material to monitor the measurements of Na<sup>+</sup>, K<sup>+</sup>, iCa<sup>2+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>. ISETROL control material is not intended for use with analyzers from other manufacturers.

**Product description**

ISETROL Electrolyte controls are available in 3 levels:

**Level 1** (red label): Low values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; high values for iCa<sup>2+</sup>.
**Level 2** (yellow label): Normal values.
**Level 3** (blue label): High values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; low values for iCa<sup>2+</sup>.

**Ingrédients**

ISETROL electrolyte control is a specially formulated non-biological solution that simulates electrolyte activity levels equivalent to 93% plasma water. ISETROL does not contain human or bovine serum proteins.

**Precautions and warnings**

For in vitro diagnostic use

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

**Storage and stability**

Unopened in the closed box up to the stated expiration date at 15-30 °C.

Do not freeze Roche Diagnostics cannot guarantee the performance of the quality controls when stored at temperatures outside the recommended range.

**FR**

**EN**

<p>ISETROL</p> Electrolyte control level 1-3	
<b>REF 03112888 180</b>	<p>10 x 1,0 mL, Niveau 1</p> 10 x 1,0 mL, Niveau 2
	<p>10 x 1,0 mL, Niveau 3</p>

**Intended use**

ISETROL, an electrolyte control is intended for use as control material to monitor the measurements of Na<sup>+</sup>, K<sup>+</sup>, iCa<sup>2+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>. ISETROL control material is not intended for use with analyzers from other manufacturers.

**Product description**

ISETROL Electrolyte controls are available in 3 levels:
**Level 1** (red label): Low values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; high values for iCa<sup>2+</sup>.
**Level 2** (yellow label): Normal values.
**Level 3** (blue label): High values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; low values for iCa<sup>2+</sup>.

**Ingrédients**

ISETROL electrolyte control is a specially formulated non-biological solution that simulates electrolyte activity levels equivalent to 93% plasma water. ISETROL does not contain human or bovine serum proteins.

**Precautions and warnings**

For in vitro diagnostic use

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

**Storage and stability**

Unopened in the closed box up to the stated expiration date at 15-30 °C.

Do not freeze Roche Diagnostics cannot guarantee the performance of the quality controls when stored at temperatures outside the recommended range.

**FR**

<p>ISETROL</p> Electrolyte control level 1-3	
<b>REF 03112888 180</b>	<p>10 x 1,0 mL, Niveau 1</p> 10 x 1,0 mL, Niveau 2
	<p>10 x 1,0 mL, Niveau 3</p>

**Intended use**

ISETROL, an electrolyte control is intended for use as control material to monitor the measurements of Na<sup>+</sup>, K<sup>+</sup>, iCa<sup>2+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>. ISETROL control material is not intended for use with analyzers from other manufacturers.

**Product description**

ISETROL Electrolyte controls are available in 3 levels:
**Level 1** (red label): Low values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; high values for iCa<sup>2+</sup>.
**Level 2** (yellow label): Normal values.
**Level 3** (blue label): High values for Na<sup>+</sup>, K<sup>+</sup>, Li<sup>+</sup> and Cl<sup>-</sup>; low values for iCa<sup>2+</sup>.

**Ingrédients**

ISETROL electrolyte control is a specially formulated non-biological solution that simulates electrolyte activity levels equivalent to 93% plasma water. ISETROL does not contain human or bovine serum proteins.

**Precautions and warnings**

For in vitro diagnostic use

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

**Storage and stability**

Unopened in the closed box up to the stated expiration date at 15-30 °C.

Do not freeze Roche Diagnostics cannot guarantee the performance of the quality controls when stored at temperatures outside the recommended range.

**FR**

<p>ISETROL</p> Electrolyte control level 1-3	
<b>REF 03112888 180</b>	<p>10 x 1,0 mL, Niveau 1</p> 10 x 1,0 mL, Niveau 2
	<p>10 x 1,0 mL, Niveau 3</p>

- Wear protective gloves, safety glasses and protective clothing.
- Prior to use gently invert the ampoule a few times to mix the contents, and to release any trapped liquid droplets. Do not shake! Avoid localized heating of ampoules.
- Protect fingers with gauze or tissue and carefully snap the neck of the ampoule.
- Aspirate liquid from the ampoule within one minute after opening, following the sampling procedure of the analyzer. Use each ampoule only once.

**Note:** Refer to the latest Instructions for use of the respective analyzer for details on the number of daily QC Measurements required for each level.
Use of heparinized sample containers for introduction of ISETROL may result in electrolyte values that are slightly different from the target values.

**Limitations**

- These controls do not contain red blood cells, and thus may not detect any system malfunctions that can cause lysis of patient red blood cells leading to erroneous results.
- These products are intended for use as quality controls, and cannot be used as calibration standards.
- The listed values only apply to the indicated lot numbers.
- These controls are used to aid in evaluation of the performance of electrolyte instrumentation, and should not be considered as a substitute for other aspects of total quality control.

**Target values and ranges**

The figures in the enclosed value sheet show the mean value and the target range of each level. The values obtained from the analyzer must fall within the defined ranges. Each defined range contains 95 % of all obtained quality control results. Refer to the Quality control section in the Instructions for use for the respective analyzer when the quality control values are outside the expected range. The target values for all levels of controls are lot-specific. Verify that the lot numbers on the value sheet match those on the ampoules.

Each laboratory must ensure that corrective measures are taken if values fall outside the limits.

Follow the applicable government regulations and local guidelines for quality control.

**Performance characteristics**

For information on performance characteristics, refer to the Instructions for use of the respective analyzer. For technical service, contact your local Roche Diagnostics organization.

**Last update:** 2013-03

FOR US CUSTOMERS ONLY: LIMITED WARRANTY Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

Significant additions or changes are indicated by a change bar in the margin.

<p>ISETROL</p> Electrolyte Control Level 1-3	
<b>REF 03112888 180</b>	<p>10 x 1,0 mL, nivel 1</p> 10 x 1,0 mL, nivel 2
	<p>10 x 1,0 mL, nivel 3</p>

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